


# D<sup>2</sup> EXTREMITIES

DESIGN AND MANUFACTURE FOR  
**HAND & FOOT SURGERY**

## HURRICANE

Stabilization screw




1.7 / 2.3 / 2.8mm



## CONICAL HEADLESS

Compression screw




**DELPHOS**  
Implants

Embracing Life™



# D<sup>2</sup> EXTREMITIES

## HURRICANE Stabilization screw




ALL SCREWS  
ARE CANNULATED

**Hurricane XXS**  
DTE-C-17XX-65H  
Diameter: 1,7mm, Cannulated: 0,9mm  
Lengths: 6-8-9-10-11-12-13-14-16-18-20-22-24-26-28-30mm

**Hurricane XS**  
DTE-C-23XX-65H  
Diameter: 2,3mm, Cannulated: 0,9mm  
Lengths: 6-8-9-10-11-12-13-14-16-18-20-22-24-26-28-30mm

**Hurricane S**  
DTE-C-28XX-65H  
Diameter: 2,8mm, Cannulated: 0,9mm / Lengths: 14-16-18-20-22-24-26-28-30mm

Four Corner Fusion



## CONICAL HEADLESS Compression screw

**Conical Headless 1,7 / DTE-C-17XX-65D**

Diameter: 1,7mm, Cannulated: 0,9mm  
Lengths: 6-7-8-9-10-11-12-13-14-15-16mm

**Conical Headless 2,3 / DTE-C-23XX-65D**

Diameter: 2,3mm, Cannulated: 0,9mm  
Lengths: 10-11-12-13-14-15-16-17-18-19-20-21-22-23-24-25-26-27-28mm

**Conical Headless 2,8 / DTE-C-28XX-65D**

Diameter: 2,8mm, Cannulated: 1,1mm  
Lengths: 10-11-12-13-14-15-16-17-18-19-20-21-22-23-24-25-26-27-28 - 30-32-34-36mm

# INSTRUMENTAL SET:

## HURRICANE Stabilization screw



**DAE-1728-865H**  
Tray for Hurricane Screws  
1.7mm, 2.3mm and 2.8mm  
*Exclusive*

## CONICAL HEADLESS Compression screw



**DAE-1728-865C**  
Tray for Conical Headless Screws  
1.7mm, 2.3mm and 2.8mm  
*Exclusive*

### Common to both system



**SK-09T06-800**  
T6 Cannulated Torx Blade



**SK-11T08-800**  
T8 Cannulated Torx Blade



**SK-C-2024-814**  
Cannulated countersink for 2.0 and 2.4mm Screws



**SK-C-1509-814**  
Cannulated Drill Ø1.5mm - Cannula 0.9 - L 90mm



**SK-C-2009-814**  
Cannulated Drill Ø2.0mm - Cannula 0.9 - L 90mm



**SK-C-2411-814**  
Cannulated Drill Ø2.4mm - Cannula 1.1 - L 90mm



**SK-0020-835**  
Tap for screws 2.0mm - L 100mm



**SK-0024-835**  
Tap for screws 2.4mm - L 100mm



**SK-0030-835**  
Tap for screws 3.0mm - L 100mm



**DA-C-0019-830**  
Cannulated screwdriver handle  
fixed axis  
(AO connections)

**DSE-0810-WIRE** Bone wire 0.8 x 100mm


**DSE-1010-WIRE** Bone wire 1.0 x 100mm



**SK-0030-800** Angled tweezer



**SK-0815-WIRE** 120mm storage tube - for 0.8mm bone wires



**SK-1015-WIRE** 120mm storage tube - for 1.0mm bone wires



**AL-0100-WIRE** Depth gauge, 100 mm, for 0.8 and 1.0 mm guide wires



**AL-0120-WIRE** Depth gauge, 100 mm, for 1.2 mm guide wires



**DSE-0000-865H**  
Sterilization Box for Hurricane Screws  
1.7mm to 2.8mm  
*Exclusive*



**DSE-0000-865C**  
Sterilization Box for Conical Headless Screws  
1.7mm to 2.8mm  
*Exclusive*

## High precision grip

Meticulous blades designed  
to make a perfect catch  
of the screw.



Embracing Life™

### Delphos Implants, S.A.

Estrada Manuel Correia Lopes,  
Parq. Empresarial Progresso, Armazém nº. 5,  
2785-126 São Domingos de Rana / Portugal

Tel.: +(351) 211 955 986

[www.delphosimplants.com](http://www.delphosimplants.com)

Our medical devices have CE Marking under TÜV Rheinland  
Notified Body, and are subjected to rigorous testing and trials.  
Delphos QMS comply with all applicable standards  
requirements of ISO 13485.