


D² EXTREMITIES

DESIGN AND MANUFACTURE FOR
HAND & FOOT SURGERY


HURRICANE
Stabilization screw



3.4 / 4.2mm



CONICAL HEADLESS
Compression screw




DELPHOS
Implants

Embracing Life™





D² EXTREMITIES


HURRICANE Stabilization screw



ALL SCREWS ARE CANNULATED

Hurricane M
DTE-C-34XX-65H  T8 
Diameter: 3,4mm, Cannulated: 1,1mm
Lengths: 14-16-18-20-22-24-26-28-30mm

Hurricane L
DTE-C-42XX-65H  T10 
Diameter: 4,2mm, Cannulated: 1,1mm
Lengths: 18-20-22-24-26-28-30-35-40-45-50mm



CONICAL HEADLESS Compression screw

Conical Headless 3,4 / DTE-C-34XX-65D
 T8 
Diameter: 3,4mm, Cannulated: 1,1mm
Lengths: 14-16-18-20-22-24-26-28-30-35-40-45-50mm

Conical Headless 4,2 / DTE-C-42XX-65D
 T8 
Diameter: 4,2mm, Cannulated: 1,1mm
Lengths: 14-16-18-20-22-24-26-28-30-35-40-45-50mm

INSTRUMENTAL SET:

HURRICANE Stabilization screw



DAE-3442-865H
Tray for Hurricane Screws
3.4mm and 4.2mm
Exclusive



SK-13T10-800
T10 Cannulated Torx Blade
Exclusive Hurricane 4.2

CONICAL HEADLESS Compression screw



DAE-3442-865C
Tray for Conical Headless Screws
3.4mm and 4.2mm
Exclusive

Common to both system



SK-11T08-800
T8 Cannulated Torx Blade



SK-C-3045-814
Cannulated countersink for 3.4 and 4.2mm Screws



SK-C-3012-814
Cannulated Drill Ø3.0mm - Cannula 1.2



SK-C-3512-814
Cannulated Drill Ø3.5mm - Cannula 1.2



SK-0035-835
Tap for screws 3.4mm - L 100mm



SK-0040-835
Tap for screws 4.2mm - L 100mm




DA-C-0019-830
Cannulated screwdriver handle - fixed axis
(AO connections)



SK-0030-800 Angled tweezer


DSE-1015-WIRE Bone wire 1.0 x 150mm




SK-1016-WIRE 160mm storage tube - for 1.0mm bone wires



AL-0120-WIRE Depth gauge, 100 mm, for 1.2 mm guide wires



DSE-0002-865H
Sterilization Box for Hurricane Screws
3.4mm to 4.2mm
Exclusive



DSE-0002-865C
Sterilization Box for Conical Headless Screws
3.4mm to 4.2mm
Exclusive

High precision grip

Meticulous blades designed to make a perfect catch of the screw.



Embracing Life™

Delphos Implants, S.A.

Estrada Manuel Correia Lopes,
Parq. Empresarial Progresso, Armazém nº. 5,
2785-126 São Domingos de Rana / Portugal

Tel.: +(351) 211 955 986

www.delphosimplants.com

Our medical devices have CE Marking under TÜV Rheinland Notified Body, and are subjected to rigorous testing and trials. Delphos QMS comply with all applicable standards requirements of ISO 13485.